

Product datasheet

Recombinant Human LHPP protein ab116175

画像数 1

製品の概要

製品名	Recombinant Human LHPP protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列

アクセッション番号 [Q9H008](#)

生物種 Human

配列 MRGSHHHHP WYASMTGGQ MGRDLYDDDD  
 KDRWGSHMAP WGKRLAGVRG VLLDISGVLY  
 DSGAGGGTAI AGSVEAVARL KRSRLKVRFC  
 TNESQKSRAE LVGQLQRLGF DISEQEV TAP  
 APAACQILKE QGLRPYLLIH DGVRSEFDQI  
 DTSNPNCVVI ADAGESFSYQ NMNNAFQVLM  
 ELEKPVLI SL GKGRYYKETS GLMLDVGPYM  
 KALEYACGIK AEVVGKPSPE FFKSALQAIG  
 VEAHQAVMIG DDIVGDVGG A QRCGMRALQV  
 RTGKFRPSDE HHPEVKADGY VDNLAEAVDL LLQHADK

分子量 34 kDa including tags

領域 1 to 270

タグ His tag N-Terminus , DDDDK tag N-Terminus

特性

Our [Abpromise guarantee](#) covers the use of **ab116175** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Mass Spectrometry
質量分析	MALDI-TOF

精製度	> 90 % SDS-PAGE. ab116175 was purified using conventional chromatography.
製品の状態	Liquid


#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 8.00 Constituents: 10% Glycerol, 0.32% Tris HCl, 0.02% DTT, 0.58% Sodium chloride
------------	--

#### 関連情報

組織特異性	Expressed in liver, kidney and moderately in brain.
配列類似性	Belongs to the HAD-like hydrolase superfamily.

#### 画像



ab116175 on a 15% SDS-PAGE (3ug)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors