

## Product datasheet

# Recombinant Human LXR beta protein ab81919

### 製品の概要

---

|        |                                    |
|--------|------------------------------------|
| 製品名    | Recombinant Human LXR beta protein |
| タンパク質長 | Full length protein                |

### 製品の詳細

---

|        |                  |
|--------|------------------|
| 由来     | Recombinant      |
| 由来     | Escherichia coli |
| アミノ酸配列 |                  |
| 生物種    | Human            |

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab81919** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生理活性 1 unit equals 1 nanogram of purified protein.

アプリケーション SDS-PAGE

精製度 > 95 % SDS-PAGE.

製品の状態 Liquid

備考 1 unit equals 1 nanogram of purified protein.

### 前処理および保存

---

保存方法および安定性 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

Preservative: None

Constituents: 20% Glycerol, 20mM Tris HCl, 100mM Potassium chloride, 1mM DTT, 0.2mM EDTA, pH 8.0

### 関連情報

---

機能 Orphan receptor. Binds preferentially to double-stranded oligonucleotide direct repeats having

the consensus half-site sequence 5'-AGGTCA-3' and 4-nt spacing (DR-4). Regulates cholesterol uptake through MYLIP-dependent ubiquitination of LDLR, VLDLR and LRP8.

**組織特異性**

Ubiquitous.

**配列類似性**

Belongs to the nuclear hormone receptor family. NR1 subfamily.  
Contains 1 nuclear receptor DNA-binding domain.

**細胞内局在**

Nucleus.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors