

Product datasheet

Recombinant Human MSRB2 protein ab101946

画像数 1

製品の概要

| | |
|--------|---------------------------------|
| 製品名 | Recombinant Human MSRB2 protein |
| タンパク質長 | Full length protein |

製品の詳細

| | |
|-----------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 由来 | Recombinant |
| 由来 | Escherichia coli |
| アミノ酸配列 | |
| アクセッション番号 | Q9Y3D2 |
| 生物種 | Human |
| 配列 | MGSSHHHHHSSGLVPRGSHMVRGQAGGGPGTGPGLGEAGSLATCELPL AKSEWQKCLTPEQFYVTREKGTPEPPFSGIYLNKEAGMYHCVCCDSPLFS SEKKYCSGTGWPSFSEAHGTSDESHTGILRRLDTSLGARTEVVCKQC EAHLGHVFPDGPNGQRFCINSVALKFKPRKH |
| 分子量 | 20 kDa including tags |
| 領域 | 21 to 182 |
| タグ | His tag N-Terminus |

特性

Our [Abpromise guarantee](#) covers the use of **ab101946** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|----------------------------------------------------------------------------------------|
| アプリケーション | SDS-PAGE Mass Spectrometry |
| 質量分析 | MALDI-TOF |
| 精製度 | > 90 % SDS-PAGE. ab101946 is purified using conventional chromatography techniques. |
| 製品の状態 | Liquid |

前処理および保存


| | |
|------------|---------------------------------------------------------------------------------------------------------|
| 保存方法および安定性 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
|------------|---------------------------------------------------------------------------------------------------------|

Preservative: None

Constituents: 10% Glycerol, 20mM Tris HCl, 1mM DTT, 0.1mM PMSF, pH 7.5

関連情報

Recombinant Human MSRB2 protein 画像



15% SDS-PAGE analysis of 3µg ab101946.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors