

## Product datasheet

# Recombinant Human MyD88 protein ab114217

画像数 1

### 製品の概要

製品名	Recombinant Human MyD88 protein
タンパク質長	Protein fragment

### 製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	<a href="#">Q99836</a>
生物種	Human
配列	RRLSLFLNVRTQVAADWTALAEEMDFEYLEIRQLETQADPTGRLLDAWQG RPGASVGRLLLELLTKLGRDDVLELGPSEEDCQKYILKQQEEAEKPLQ
分子量	37 kDa including tags
領域	31 to 130

### 特性

Our [Abpromise guarantee](#) covers the use of **ab114217** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA SDS-PAGE Western blot
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml


### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

## 関連情報

機能	Adapter protein involved in the Toll-like receptor and IL-1 receptor signaling pathway in the innate immune response. Acts via IRAK1, IRAK2, IRF7 and TRAF6, leading to NF-kappa-B activation, cytokine secretion and the inflammatory response. Increases IL-8 transcription. Involved in IL-18-mediated signaling pathway.
組織特異性	Ubiquitous.
関連疾患	Defects in MYD88 are the cause of MYD88 deficiency (MYD88D) [MIM:612260]; also known as recurrent pyogenic bacterial infections due to MYD88 deficiency. Patients suffer from autosomal recessive, life-threatening, often recurrent pyogenic bacterial infections, including invasive pneumococcal disease, and die between 1 and 11 months of age. Surviving patients are otherwise healthy, with normal resistance to other microbes, and their clinical status improved with age.
配列類似性	Contains 1 death domain. Contains 1 TIR domain.
ドメイン	The intermediate domain (ID) is required for the phosphorylation and activation of IRAK.
細胞内局在	Cytoplasm.

## 画像



ab114217 analysed on a 12.5% SDS-PAGE gel stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors