

## Product datasheet

# Recombinant Human PAFAH protein ab87477

### 画像数 1

#### 製品の概要

---

|        |                                 |
|--------|---------------------------------|
| 製品名    | Recombinant Human PAFAH protein |
| タンパク質長 | Full length protein             |

#### 製品の詳細

---

|        |                  |
|--------|------------------|
| 由来     | Recombinant      |
| 由来     | Escherichia coli |
| アミノ酸配列 |                  |
| 生物種    | Human            |

#### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab87477** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

|          |                  |
|----------|------------------|
| アプリケーション | SDS-PAGE         |
| 精製度      | > 95 % SDS-PAGE. |
| 製品の状態    | Liquid           |

#### 前処理および保存

---

|            |  |
|------------|--|
| 保存方法および安定性 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.<br>Preservative: None<br>Constituents: 50% Glycerol, PBS |
|------------|--|

#### 関連情報

---

|       |   |
|-------|---|
| 機能    | Modulates the action of platelet-activating factor (PAF) by hydrolyzing the sn-2 ester bond to yield the biologically inactive lyso-PAF. Has a specificity for substrates with a short residue at the sn-2 position. It is inactive against long-chain phospholipids. |
| 組織特異性 | Plasma.   |

## 関連疾患

Platelet-activating factor acetylhydrolase deficiency

Asthma

Atopic hypersensitivity


## 配列類似性

Belongs to the AB hydrolase superfamily. Lipase family.

## 細胞内局在

Secreted, extracellular space.

## 画像



ab87477 on SDS-PAGE.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors