

Product datasheet

Recombinant Human PC6 protein ab114944

画像数 1

製品の概要

製品名	Recombinant Human PC6 protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	Q92824
生物種	Human
配列	TEGYFMEDGRCVQSCSISYYFDHSSSENGYKSCKKCDISCLTCNGPGFKNC TSCPSGYLLDLGMCQMGAICKDATEESWAEGGFCMLVKNNLCQRKVLQQ LCCKTCTFQG
分子量	38 kDa including tags
領域	804 to 913

特性

Our [Abpromise guarantee](#) covers the use of **ab114944** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA SDS-PAGE Western blot
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml. Best used within three months from the date of receipt.

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------	-------------------------------------------------------------------------------------------


pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報

機能	Likely to represent a widespread endoprotease activity within the constitutive and regulated secretory pathway. Capable of cleavage at the RX(K/R)R consensus motif.
組織特異性	Expressed in T-lymphocytes.
配列類似性	Belongs to the peptidase S8 family. Contains 1 homo B/P domain. Contains 1 PLAC domain.
ドメイン	The propeptide domain acts as an intramolecular chaperone assisting the folding of the zymogen within the endoplasmic reticulum.
細胞内局在	Secreted.

画像



ab114944 analyzed on a 12.5% SDS-PAGE gel stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors