

## Product datasheet

# Recombinant Human PCNP protein ab113147

### 画像数 1

#### 製品の概要

製品名	Recombinant Human PCNP protein
タンパク質長	Full length protein

#### 製品の詳細

由来	Recombinant
由来	Escherichia coli

#### アミノ酸配列

アクセッション番号 [Q8WW12](#)

生物種 Human

配列 **MGSSHHHHHSSGLVPRGSHMGSHMADGKAGDEKPEKSQRAGAAGGPEEE**  
**AEKPVKTKTVSSNGGESSRSAEKRSAEAAAADLPTKPTKISKFGFAIG**  
**SQTTKKASAIKLGSSPKETVPTLAPKTLVAAAFNEDEDESEPEEMPP**  
**EAKMRMKNIGRDTPTSAGPNSFNKGKHGFSNQKLWERNIKSHLGNVHDQ**  
**DN**

分子量	22 kDa including tags
領域	1 to 178
タグ	His tag N-Terminus

#### 特性

Our [Abpromise guarantee](#) covers the use of **ab113147** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Mass Spectrometry SDS-PAGE
質量分析	MALDI-TOF
精製度	> 95 % SDS-PAGE. ab113147 was purified using conventional chromatography.
製品の状態	Liquid

## 保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.32% Tris HCl, 10% Glycerol, 0.02% DTT

## 関連情報

### 機能

May be involved in cell cycle regulation.


### 翻訳後修飾

Ubiquitinated; mediated by UHRF2 and leading to its subsequent proteasomal degradation.

### 細胞内局在

Nucleus.

## 画像



15% SDS-PAGE analysis of ab113147 (3µg)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors