

Product datasheet

Recombinant Human PSMG4 protein ab117216

画像数 1

製品の概要

製品名	Recombinant Human PSMG4 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli
アミノ酸配列	
アクセッション番号	Q5JS54
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMEGLVVAAGGDVSLHNFSARLWEQLVHFHV MRLTDSLFLWVGATPHLRNLAVAMCSRYDSIPVSTSLLGDTSDTTSTGLA QRLARKTNKQVFVSYNLQNTDSNFALLVENRIKEEMAEAFPEKF
分子量	16 kDa including tags
領域	1 to 123
タグ	His tag N-Terminus

特性

Our [Abpromise guarantee](#) covers the use of **ab117216** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Mass Spectrometry SDS-PAGE
質量分析	MALDI-TOF
精製度	> 95 % SDS-PAGE. ab117216 was purified using conventional chromatography techniques.
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00


Constituents: 0.32% Tris HCl, 0.02% DTT, 10% Glycerol, 0.58% Sodium chloride

関連情報

機能 Chaperone protein which promotes assembly of the 20S proteasome.

配列類似性 Belongs to the PSMG4 family.

画像



15% SDS-PAGE analysis of ab117216 (3µg)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors