

Product datasheet

Recombinant Human Pallidin protein ab124579

画像数 1

製品の概要

製品名	Recombinant Human Pallidin protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列

アクセッション番号 [Q9UL45](#)

生物種 Human

配列
 MGSSHHHHHH SGLVPRGSH MSVPGPSSPD
 GALTRPPYCL EAGEPTPLS DTSPDEGLIE
 DLTIEDKAVE QLAEGLLSHY LPDLQRSKQA
 LQELTQNQVV LLDTLEQEIS KFKECHSMLD
 INALFAEAKH YHAKLVNIRK EMLMLHEKTS
 KLKKRALKLQ QKRQKEELER EQQREKEFER
 EKQLTARPAK RM

分子量 22 kDa including tags

領域 1 to 172

タグ His tag N-Terminus

特性

Our [Abpromise guarantee](#) covers the use of **ab124579** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
精製度	> 90 % SDS-PAGE. ab124579 is purified using conventional chromatography techniques.
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.03% DTT, 0.32% Tris HCl, 10% Glycerol, 0.58% Sodium chloride

関連情報

機能

Involved in the development of lysosome-related organelles, such as melanosomes and platelet-dense granules. May play a role in intracellular vesicle trafficking, particularly in the vesicle-docking and fusion process.

組織特異性

Widely expressed.


翻訳後修飾

Phosphorylated.

細胞内局在

Cytoplasm. Endomembrane system. It can exist as a soluble protein as well as a peripheral membrane protein.

画像



15% SDS-PAGE showing ab124579 (3µg).

SDS-PAGE - Pallidin protein (His tag) (ab124579)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors