

Product datasheet

Recombinant Human RABGGTA protein ab116789

画像数 1

製品の概要

製品名	Recombinant Human RABGGTA protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Wheat germ

アミノ酸配列

アクセッション番号 [Q92696](#)

生物種 Human

配列
 MHGRLKVKTSEEQAEAKRLEREQKLLKLYQSATQAVFQKRQAGELDESIVLE
 LTSQILGANPDFATLWNCRRREVLQQLETQKSPEELAALVKAELGFLESCL
 RVNPKSYGTWHHRCWLLGRLPEPNWTRELELCARFLEVDERNFHCWDYRR
 FVATQAAVPPAEELAFDLSLITRNFSNYSSWHYRSCLLPQLHPQPDSGPQ
 GRLPEDVLLKELELVQNAFFDTPNDQSAWFYHRWLLGRADPQDALRCLHV
 SRDEACLTVSFSRPLLVGSRMEILLMVDDSPILIVEWRTPDGRNRPSHVW
 LCDLPAASLNDQLPQHTFRVIWTAGDVQKECVLLKGRQEGWCRDSTTDEQ
 LFRCELSVEKSTVLQSELESCKELQELEPENKWCLLTIIILLMRALDPLLY
 EKETLQYFQTLKAVDPMRATYLDDLRSKFLENSVLKMEYAEVRVHLAH
 KDLTVLCHLEQLLLVTHLDLSHNRLRTLPPALAALRCLEVLQASDNAIES
 LDGVTNLPRLQELLLCNRLQQPAVLQPLASCPRLVLLNLQGNPLCQAVG
 ILEQLAELLPSVSSVLT

分子量 88 kDa including tags

領域 1 to 567

特性

Our [Abpromise guarantee](#) covers the use of **ab116789** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE
	ELISpot

製品の状態

Liquid

備考

Protein concentration is above or equal to 0.05 mg/ml.
Best used within three months from the date of receipt.

前処理および保存**保存方法および安定性**

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報**機能**


Catalyzes the transfer of a geranyl-geranyl moiety from geranyl-geranyl pyrophosphate to both cysteines in Rab proteins with an -XXCC, -XCXC and -CCXX C-terminal, such as RAB1A, RAB3A and RAB5A respectively.

配列類似性

Belongs to the protein prenyltransferase subunit alpha family.

Contains 5 LRR (leucine-rich) repeats.

Contains 6 PFTA repeats.

画像

12.5% SDS-PAGE Stained with Coomassie Blue.

SDS-PAGE - RABGGTA protein (ab116789)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors