

Product datasheet

Recombinant Human RABGGTA protein ab116789

画像数 1

製品の概要

製品名	Recombinant Human RABGGTA protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	<a href="#">Q92696</a>
生物種	Human
配列	MHGRLKVKTSEEQAEAKRLEREQKLLKYQSATQAVFQKRQAGELDESVLE LTSQILGANPDFATLWNCRRREVLQQLTQKSPEELAALVKAELGFLESCL RVNPKSYGTWHHRCWLLGRLPEPNWTRELELCARFLEVDERNFHCWDYRR FVATQAAVPPAEELAFDLSLITRNFSNYSSWHYRSCLLPQLHPQPDSGPQ GRLPEDVLLKELELVQNAFFTPNDQSAWFYHRWLLGRADPQDALRCLHV SRDEACLTVSFSRPLLVGSRMEILLMVDDSPDIVERTPDGRNRPSHW LCDLPAASLNDQLPQHTFRVIWTAGDVQKECVLLKGRQEGWCRDSTTDEQ LFRCELSVEKSTVLQSELESCKELQELEPENKWCLLTIILLMRALDPLLY EKETLQYFQTLKAVDPMRATYLDDLRSKFLENSVLKMEYAEVRVLAHLA KDLTVLCHLEQLLLVTHLDLSHNRLRTLPPALAALRCLEVLQASDNAIES LDGVTNLPRLQELLLCNRLQQPAVLQPLASCPRLVLLNLQGNPLCQAVG ILEQLAELLPSVSSVLT
分子量	88 kDa including tags
領域	1 to 567

特性

Our [Abpromise guarantee](#) covers the use of **ab116789** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE ELISpot
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml.

## 前処理および保存

### 保存方法および安定性

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

## 関連情報

### 機能

Catalyzes the transfer of a geranyl-geranyl moiety from geranyl-geranyl pyrophosphate to both cysteines in Rab proteins with an -XXCC, -XCXC and -CCXX C-terminal, such as RAB1A, RAB3A and RAB5A respectively.


### 配列類似性

Belongs to the protein prenyltransferase subunit alpha family.

Contains 5 LRR (leucine-rich) repeats.

Contains 6 PFTA repeats.

## Recombinant Human RABGGTA protein 画像



12.5% SDS-PAGE Stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors