

Product datasheet

Recombinant Human RABGGTA protein ab116789

画像数 1

製品の概要

製品名	Recombinant Human RABGGTA protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Wheat germ

アミノ酸配列

アクセッション番号 [Q92696](#)

生物種 Human

配列  
 MHGRLKVKTSEEQAEAKRLEREQKLLKLYQSATQAVFQKRQAGELDESIVLE  
 LTSQILGANPDFATLWNCRRREVLQQLETQKSPEELAALVKAELGFLESCL  
 RVNPKSYGTWHHRCWLLGRLPEPNWTRELELCARFLEVDERNFHCWDYRR  
 FVATQAAVPPAEELAFDLSLITRNFSNYSSWHYRSCLLPQLHPQPDSGPQ  
 GRLPEDVLLKELELVQNAFFDTPNDQSAWFYHRWLLGRADPQDALRCLHV  
 SRDEACLTVSFSRPLLVGSRMEILLMVDDSPDIVERTPDGRNRPSHW  
 LCDLPAASLNDQLPQHTFRVIWTAGDVQKECVLLKGRQEGWCRDSTTDEQ  
 LFRCELSVEKSTVLQSELESCKELQELEPENKWCLLTIILLMRALDPLLY  
 EKETLQYFQTLKAVDPMRATYLDDLRSKFLENSVLKMEYAEVRVHLAH  
 KDLTVLCHLEQLLLVTHLDLSHNRLRTLPPALAALRCLEVLQASDNAIES  
 LDGVTNLPRLQELLLCNRLQQPAVLQPLASCPRLVLLNLQGNPLCQAVG  
 ILEQLAELLPSVSSVLT

分子量 88 kDa including tags

領域 1 to 567

特性

Our [Abpromise guarantee](#) covers the use of **ab116789** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE
	ELISpot

**製品の状態**

Liquid

**備考**

Protein concentration is above or equal to 0.05 mg/ml.  
Best used within three months from the date of receipt.

**前処理および保存****保存方法および安定性**

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

**関連情報****機能**


Catalyzes the transfer of a geranyl-geranyl moiety from geranyl-geranyl pyrophosphate to both cysteines in Rab proteins with an -XXCC, -XCXC and -CCXX C-terminal, such as RAB1A, RAB3A and RAB5A respectively.

**配列類似性**

Belongs to the protein prenyltransferase subunit alpha family.

Contains 5 LRR (leucine-rich) repeats.

Contains 6 PFTA repeats.

**画像**

12.5% SDS-PAGE Stained with Coomassie Blue.

SDS-PAGE - RABGGTA protein (ab116789)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors