

Product datasheet

Recombinant Human SNAPC4 protein ab116993

画像数 1

製品の概要

| | |
|--------|----------------------------------|
| 製品名 | Recombinant Human SNAPC4 protein |
| タンパク質長 | Protein fragment |

製品の詳細

| | |
|----|-------------|
| 由来 | Recombinant |
| 由来 | Wheat germ |

アミノ酸配列

| | |
|-----------|--|
| アクセッション番号 | Q5SXM2 |
| 生物種 | Human |
| 配列 | PADPPISEEERWGEASNDEDDPKDKTLPEDPETCLQLNMVYQEVIQEKLA EANLLLAQNREQEELMRDLAGSKGKVKDGGKSLPPSTYMGHFMPYFKD KVTGVGPPAN |
| 分子量 | 38 kDa including tags |
| 領域 | 53 to 162 |

特性

Our [Abpromise guarantee](#) covers the use of **ab116993** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|---|
| アプリケーション | SDS-PAGE ELISA Western blot |
| 製品の状態 | Liquid |
| 備考 | Protein concentration is above or equal to 0.05 mg/ml. This protein is best used within three months from the date of receipt. |

前処理および保存


| | |
|------------|---|
| 保存方法および安定性 | Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. |
|------------|---|

pH: 8.00
Constituents: 0.79% Tris HCl, 0.3% Glutathione

関連情報

| | |
|-------|---|
| 機能 | Part of the SNAPc complex required for the transcription of both RNA polymerase II and III small-nuclear RNA genes. Binds to the proximal sequence element (PSE), a non-TATA-box basal promoter element common to these 2 types of genes. Recruits TBP and BRF2 to the U6 snRNA TATA box. |
| 配列類似性 | Contains 3 HTH myb-type DNA-binding domains. Contains 2 Myb-like domains. |
| 細胞内局在 | Nucleus. |

画像



12.5% SDS-PAGE showing ab116993 at approximately 37.84kDa.
Stained with Coomassie Blue.

SDS-PAGE - SNAPC4 protein (Human) (ab116993)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors