

Product datasheet

Recombinant Human SNAPC4 protein ab116993

画像数 1

製品の概要

製品名	Recombinant Human SNAPC4 protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	Q5SXM2
生物種	Human
配列	PADPPPISEEERWGEASNDEDDPKDKTLPEDPETCLQLNMVYQEVIQEKLA EANLLLAQNREQQEELMRDLAGSKGTKVKDGLKSLPPSTYMGHFMKPYFKD KVTGVGPPAN
分子量	38 kDa including tags
領域	53 to 162

特性

Our [Abpromise guarantee](#) covers the use of **ab116993** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE ELISA Western blot
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml. This protein is best used within three months from the date of receipt.


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.79% Tris HCl, 0.3% Glutathione
------------	---

関連情報

機能	Part of the SNAPc complex required for the transcription of both RNA polymerase II and III small-nuclear RNA genes. Binds to the proximal sequence element (PSE), a non-TATA-box basal promoter element common to these 2 types of genes. Recruits TBP and BRF2 to the U6 snRNA TATA box.
配列類似性	Contains 3 HTH myb-type DNA-binding domains. Contains 2 Myb-like domains.
細胞内局在	Nucleus.

Recombinant Human SNAPC4 protein 画像



12.5% SDS-PAGE showing ab116993 at approximately 37.84kDa.
Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors