

Product datasheet

Recombinant Human SP2 transcription factor protein ab117031

画像数 1

製品の概要

製品名	Recombinant Human SP2 transcription factor protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Wheat germ

アミノ酸配列

アクセッション番号 [Q02086](#)

生物種 Human

配列
 MAATAAVSPSDYLQPAASTTQDSQPSPLALLAATCSKIGPPAVEAAVTPP
 APPQPTPRKLVPIKPAPLPLSPGKNSFGILSSKGNILQIQGSQLSASYPG
 GQLVFAIQNPTMINKGTRSNANIQYQAVPQIQASNSQTIQVQPNLTNIIQ
 IIPGTNQAIITPSPSSHKPVPIKPAPIQKSSTTTTPVQSGANVVKLTGGG
 GNVTLTLPVNNLVNASDTGAPTQLLTASCQTGMLNSTRMFLFLAFINVL

分子量 53 kDa including tags

領域 1 to 249

特性

Our [Abpromise guarantee](#) covers the use of **ab117031** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	SDS-PAGE
	Western blot

製品の状態 Liquid

備考 Protein concentration is above or equal to 0.05 mg/ml.
 Best use within three months from the date of receipt of this protein.

前処理および保存

保存方法および安定性

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報

機能

Binds to GC box promoters elements and selectively activates mRNA synthesis from genes that contain functional recognition sites.


配列類似性

Belongs to the Sp1 C2H2-type zinc-finger protein family.
Contains 3 C2H2-type zinc fingers.

細胞内局在

Nucleus.

画像



12.5% SDS-PAGE stained with Coomassie

Blue showing ab117031 at approximately
53.46 kDa.

SDS-PAGE - SP2 transcription factor protein
(Human) (ab117031)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors