

Product datasheet

Recombinant Human Slit3 protein ab116753

画像数 1

製品の概要

製品名	Recombinant Human Slit3 protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	<a href="#">O75094</a>
生物種	Human
配列	PCLGHRCHHGKCVATGTSYMCKCAEGYGGDLCDNKNSANACSAFKCHHG QCHISDQGEPEYCLCQPGFSGEHCQQENPCLGQVVREVIRRQKGYASCATA
分子量	37 kDa including tags
領域	1371 to 1470

特性

Our [Abpromise guarantee](#) covers the use of **ab116753** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA SDS-PAGE Western blot
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml. Best use within three months from the date of receipt of this protein.

前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------	---

Constituents: 0.79% Tris HCl, 0.3% Glutathione

## 関連情報

機能	May act as molecular guidance cue in cellular migration, and function may be mediated by interaction with roundabout homolog receptors.
組織特異性	Predominantly expressed in thyroid.
配列類似性	Contains 1 CTCK (C-terminal cystine knot-like) domain. Contains 9 EGF-like domains. Contains 1 laminin G-like domain. Contains 20 LRR (leucine-rich) repeats. Contains 4 LRRCT domains. Contains 4 LRRNT domains.
細胞内局在	Secreted.

## 画像



12.5% SDS-PAGE analysis of ab116753 at approximately 36.63kDa, stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors