

Product datasheet

Recombinant Human Slit3 protein ab116753

画像数 1

製品の概要

製品名	Recombinant Human Slit3 protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Wheat germ
アミノ酸配列	
アクセッション番号	O75094
生物種	Human
配列	PCLGHRCHHGKCVATGTSYMCKCAEGYGGDLCDNKNDANACSAFKCHHG QCHISDQGEPEYCLCQPGFSGEHCQENPCLGQVVREVI RRQKGYASCATA
分子量	37 kDa including tags
領域	1371 to 1470

特性

Our [Abpromise guarantee](#) covers the use of **ab116753** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA SDS-PAGE Western blot
製品の状態	Liquid
備考	Protein concentration is above or equal to 0.05 mg/ml. Best use within three months from the date of receipt of this protein.

前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------	---

Constituents: 0.79% Tris HCl, 0.3% Glutathione

関連情報

機能	May act as molecular guidance cue in cellular migration, and function may be mediated by interaction with roundabout homolog receptors.
組織特異性	Predominantly expressed in thyroid.
配列類似性	Contains 1 CTCK (C-terminal cystine knot-like) domain. Contains 9 EGF-like domains. Contains 1 laminin G-like domain. Contains 20 LRR (leucine-rich) repeats. Contains 4 LRRCT domains. Contains 4 LRRNT domains.
細胞内局在	Secreted.

画像



12.5% SDS-PAGE analysis of ab116753 at approximately 36.63kDa, stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors