

Product datasheet

Recombinant Human TNFSF18 protein ab113132

画像数 1

製品の概要

製品名	Recombinant Human TNFSF18 protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Escherichia coli
アミノ酸配列	
アクセッション番号	Q9UNG2
生物種	Human
配列	MQLETAKEPCMAKFGPLPSKWQMASSEPPCVNKVSDWKLEILQNGLYLIY GQVAPNANYNDVAPFEVRLYKNKDMIQTLTNKSKIQNVGGTYELHVGD TI DLIFNSEHQVLKNNNTYWGIILLANPQFIS
分子量	15 kDa
領域	72 to 199

特性

Our [Abpromise guarantee](#) covers the use of **ab113132** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
精製度	> 90 % SDS-PAGE. ab113132 was purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer.
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------	--


pH: 3.50

Constituents: 0.26% Sodium citrate, 10% Glycerol, 0.02% DTT

関連情報

機能	Cytokine that binds to TNFRSF18/AITR/GITR. Important for interactions between activated T-lymphocytes and endothelial cells and may modulate T-lymphocyte survival in peripheral tissues.
組織特異性	Expressed at high levels in the small intestine, ovary, testis, kidney and endothelial cells.
配列類似性	Belongs to the tumor necrosis factor family.
細胞内局在	Membrane.

Recombinant Human TNFSF18 protein 画像



15% SDS-PAGE analysis of 3 µg ab113132.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors