

## Product datasheet

# Recombinant Human TNFSF18 protein ab113132

### 画像数 1

#### 製品の概要

製品名	Recombinant Human TNFSF18 protein
タンパク質長	Protein fragment

#### 製品の詳細

由来	Recombinant
由来	Escherichia coli

#### アミノ酸配列

アクセッション番号	<a href="#">Q9UNG2</a>
生物種	Human
配列	MQLETAKEPCMAKFGPLPSKWQMASSEPPCVNKVSDWKLEILQNGLYLIY GQVAPNANYNDVAPFEVRLYKNKDMIQTLTNKSKIQNVGGTYELHVG D TI DLIFNSEHQVLKNNNTYWGIILLANPQFIS
分子量	15 kDa
領域	72 to 199

#### 特性

Our [Abpromise guarantee](#) covers the use of **ab113132** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
精製度	> 90 % SDS-PAGE. ab113132 was purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer.
製品の状態	Liquid

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -
------------	--------------------------------------------------------------------------------------------------

80°C. Avoid freeze / thaw cycle.


pH: 3.50

Constituents: 0.26% Sodium citrate, 10% Glycerol, 0.02% DTT

## 関連情報

機能	Cytokine that binds to TNFRSF18/AITR/GITR. Important for interactions between activated T-lymphocytes and endothelial cells and may modulate T-lymphocyte survival in peripheral tissues.
組織特異性	Expressed at high levels in the small intestine, ovary, testis, kidney and endothelial cells.
配列類似性	Belongs to the tumor necrosis factor family.
細胞内局在	Membrane.

## 画像



15% SDS-PAGE analysis of 3 µg ab113132.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors