

Product datasheet

Recombinant Human VAPA protein ab93462

画像数 1

製品の概要

製品名	Recombinant Human VAPA protein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Escherichia coli

アミノ酸配列

生物種	Human
配列	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHTMAS ASGAMAKHEQ ILVLDPPTDL KFKGPFTDVV TTNLKLRNPS DRKVCFKVKT TAPRRYCVRP NSGIIDPGST VTVSVMLQPF DYDPNEKSKH KFMVQTI FAP PNTSDMEAVW KEAKPDELMD SKLRVCFEMP NENDKLN DME PSKAVPLNAS KQDGMPKPH SVSLNDTETR KLMEECKRLQ GEMMKLSEEN RHLRDEGLRL RKVAHSDKPG STSTASFRDN VTSP

領域	1 to 227
----	----------

特性

Our [Abpromise guarantee](#) covers the use of **ab93462** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
----------	----------

製品の状態	Liquid
-------	--------

前処理および保存


保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Preservative: None
------------	---

Constituents: 10% Glycerol, 20mM Tris HCl, 1mM DTT, pH 8.0

関連情報

機能	May play a role in vesicle trafficking.
組織特異性	Ubiquitous.
配列類似性	Belongs to the VAMP-associated protein (VAP) (TC 9.B.17) family. Contains 1 MSP domain.
細胞内局在	Endoplasmic reticulum membrane. Present in the plasma membrane and in intracellular vesicles, together with SNARE proteins. May also associate with the cytoskeleton. Colocalizes with OCLN at the tight junction in polarized epithelial cells.

画像



15% SDS-PAGE analysis of ab93462 (3µg).

SDS-PAGE - VAPA protein (His tag) (ab93462)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors