

## Product datasheet

# Recombinant Mouse IL28A protein ab50172

### 製品の概要

---

製品名	Recombinant Mouse IL28A protein
タンパク質長	Protein fragment

### 製品の詳細

---

由来	Recombinant
由来	Escherichia coli

### アミノ酸配列

生物種	Mouse
配列	MDPVPRATRL PVEAKDCHIA QFKSLSPKEL QAFKKAKDAI EKRLLEKDMR CSSHLISRAW DLKQLQVQER PKALQAEVAL TLKVWENMTD SALATILGQP LHTLSHIHSQ LQTCTQLQAT AEPKPPSRRL SRWLHRLQEA QSKETPGCLE DSVTSNLFRL LTRDLKCVAS GDQCV

分子量	20 kDa
領域	19 to 193

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab50172** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
エンドトキシン・レベル	< 0.100 Eu/μg
精製度	> 95 % SDS-PAGE. ab50172 purity was assessed also by HPLC Endotoxin level is less than 0.1 ng per μg (1EU/μg).
製品の状態	Lyophilised

### 前処理および保存

---

<b>保存方法および安定性</b>	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.
<b>再構成</b>	Reconstitute to 1mg/ml in distilled water.
<b>関連情報</b>	
<b>機能</b>	Cytokine with immunomodulatory activity. Up-regulates MHC class I antigen expression. Displays potent antiviral activity. Also displays antitumor activity. Ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IL28RA. The ligand/receptor complex seems to signal through the Jak-STAT pathway.
<b>配列類似性</b>	Belongs to the IL-28/IL-29 family.
<b>細胞内局在</b>	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors