

Product datasheet

Recombinant Rat Resistin protein ab52152

製品の概要

製品名	Recombinant Rat Resistin protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli
アミノ酸配列	
生物種	Rat

特性

Our [Abpromise guarantee](#) covers the use of **ab52152** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot ELISA
精製度	> 95 % SDS-PAGE. >95% by SDS-PAGE.
製品の状態	Lyophilised
備考	The reconstituted solution can be diluted into aqueous buffers and stored at 4°C for 1 week or –20°C for future use.

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. 20mM Tris, pH 8 Endotoxin <0.1 ng/μg (IEU/μg).
再構成	Reconstitute in water to a concentration of 0.1-1.0 mg/ml.

関連情報

機能	Hormone that seems to suppress insulin ability to stimulate glucose uptake into adipose cells. Potentially links obesity to diabetes.
組織特異性	Expressed only in fatty tissues.
配列類似性	Belongs to the resistin/FIZZ family.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors