

Product datasheet

Recombinant SARS spike glycoprotein ab49044

製品の概要

製品名	Recombinant SARS spike glycoprotein
タンパク質長	Protein fragment

製品の詳細

由来	Recombinant
由来	Escherichia coli

特性

Our [Abpromise guarantee](#) covers the use of **ab49044** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	Western blot

製品の状態	Liquid
-------	--------

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
	Preservative: 0.01% Sodium Azide
	Constituents: 50% Glycerol, 1.5M Urea, 25mM Tris HCl, 1mM EDTA

関連情報

関連性	A novel coronavirus has been identified as the causative agent of SARS (Severe Acute Respiratory Syndrome). Coronaviruses are a major cause of upper respiratory diseases in humans. The genomes of these viruses are positive stranded RNA approximately 27 to 31kb in length. SARS infection can be mediated by the binding of the viral spike protein, a glycosylated 139 kDa protein and the major surface antigen of the virus, to the angiotensin converting enzyme 2 (ACE2) on target cells. This binding can be blocked by a soluble form of ACE2.
-----	--

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors