

Product datasheet

Recombinant *Treponema pallidum* p47 protein (Beta-Galactosidase) ab51323

製品の概要

製品名	Recombinant <i>Treponema pallidum</i> p47 protein (Beta-Galactosidase)
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	<i>Escherichia coli</i>
アミノ酸配列	
配列の追加情報	Fused with β -galactosidase at the N-terminus.

特性

Our [Abpromise guarantee](#) covers the use of **ab51323** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生理活性	Reacts strongly with human <i>T. pallidum</i> positive serum.
アプリケーション	SDS-PAGE Flow Cytometry Western blot ELISA
精製度	> 95 % SDS-PAGE. Purity of proteins is evaluated by SDS-PAGE, by measuring optical density at 280 nm and by method of Bradford et al.
製品の状態	Liquid
備考	Reacts strongly with human <i>T. pallidum</i> positive serum.

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Preservative: None
------------	---

Constituents: 8M Urea, 20mM Tris HCl, 10mM Beta mercaptoethanol, pH 8

関連情報

関連性

Treponema pallidum is the causative agent of syphilis. It is a gram negative spirochete, a helical to sinusoidal bacterium with outer and cytoplasmic membranes, a thin peptidoglycan layer, and periplasmic flagella.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors