

## Product datasheet

# Recombinant human C5a protein ab61918

### 製品の概要

---

製品名	Recombinant human C5a protein
タンパク質長	Full length protein

### 製品の詳細

---

由来	Recombinant
由来	Escherichia coli

### アミノ酸配列

生物種	Human
配列	TLQKKIEEIA AKYKHSVVKK CCYDGACVNN DETCEQRAAR ISLGPRCIKA FTECCVVASQ LRANISHKDM QLGR

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab61918** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
精製度	> 95 % SDS-PAGE.
製品の状態	Lyophilised

### 前処理および保存

---

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Preservative: None Constituents: 2.5mg/ml BSA This product is an active protein and may elicit a biological response in vivo, handle with caution.
再構成	Reconstitute to a concentration of 1.0 mg/ml in ice-cold water containing 2.5 mg/ml BSA.

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors