

## Product datasheet

# Recombinant human G-CSF protein ab9692

### 製品の概要

---

製品名	Recombinant human G-CSF protein
タンパク質長	Full length protein

### 製品の詳細

---

由来	Recombinant
由来	Escherichia coli

### アミノ酸配列

生物種	Human
配列	Human G-CSF is a 18.7 kDa protein containing 174 amino acid residues: TPLGPASSLP QSFLLKCLEQ VRKIQGDGAA LQEKLCATYK LCHPEELVLL GHSLG IPWAPLSSCPSQALQ LAGCLSQLHS GLFLYQGLLQ ALEGISPELG P TLDTLQLDV ADFATTIWQQMEELGMAPAL QPTQGAMPAF ASAFQRRRA GG VLVASHLQSF LEVSYRVLRH LAQP

### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab9692** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生理活性	The ED <sub>50</sub> , determined by the dose-dependent stimulation of the proliferation of murine NSF-60 cells, is < 0.1 ng/ml, corresponding to a specific activity of > 1 x 10 <sup>7</sup> units/mg.
アプリケーション	Functional Studies SDS-PAGE
エンドトキシン・レベル	< 0.100 Eu/μg
精製度	> 95 % SDS-PAGE. Sterile filtered Greater than 98% pure by SDS-PAGE and HPLC analyses.
製品の状態	Lyophilised
備考	The ED <sub>50</sub> as determined by the dose-dependent stimulation of the proliferation of murine NFS-60 cells is ≤ 0.1 ng/ml, corresponding to a specific activity of ≥ 1 x 10 <sup>7</sup> units/mg.

## 前処理および保存

---

### 保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

n/a

This product is an active protein and may elicit a biological response in vivo, handle with caution.

## 関連情報

---

### 機能

Granulocyte/macrophage colony-stimulating factors are cytokines that act in hematopoiesis by controlling the production, differentiation, and function of 2 related white cell populations of the blood, the granulocytes and the monocytes-macrophages. This CSF induces granulocytes.

### 配列類似性

Belongs to the IL-6 superfamily.

### 翻訳後修飾

O-glycan consists of Gal-GalNAc disaccharide which can be modified with up to two sialic acid residues (done in recombinantly expressed G-CSF from CHO cells).

### 細胞内局在

Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors