

Product datasheet

Recombinant human RAB14 protein ab90105

製品の概要

製品名	Recombinant human RAB14 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Escherichia coli
アミノ酸配列	
生物種	Human

特性

Our [Abpromise guarantee](#) covers the use of **ab90105** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生理活性	Activity: 100 pmol of protein can bind > 80 pmol of GDP.
アプリケーション	Functional Studies SDS-PAGE
精製度	> 90 % SDS-PAGE.
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 50mM Tris HCl, 2mM Beta mercaptoethanol, 100mM Sodium chloride, 10mM Magnesium chloride, 0.01mM GDP, pH 8.0 This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

関連情報

機能	May be involved in vesicular trafficking and neurotransmitter release.
配列類似性	Belongs to the small GTPase superfamily. Rab family.
細胞内局在	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors