

Product datasheet

Anti-S100A4 antibody, prediluted ab27427

1 References [画像数 2](#)

製品の概要

製品名	Anti-S100A4 antibody, prediluted
製品の詳細	Rabbit polyclonal to S100A4, prediluted
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB
種交差性	交差種: Human
免疫原	Recombinant protein encoding full length human S100A4.
ポジティブ・コントロール	Tonsil and melanoma.

製品の特性

製品の状態	Prediluted
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.1% Sodium Azide Constituents: PBS, BSA
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab27427** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/1. Use neat for 10 min at room temperature. Perform heat mediated antigen retrieval via the pressure cooker method before commencing with IHC staining protocol, by boiling tissue sections in 10mM citrate buffer, pH 6.0 for 10 min followed by cooling at room temperature for 20 min.
WB		Use at an assay dependent concentration. Predicted molecular weight: 12 kDa.


ターゲット情報

組織特異性	Ubiquitously expressed.
配列類似性	Belongs to the S-100 family. Contains 2 EF-hand domains.

画像

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-S100A4 antibody, prediluted (ab27427)

Ab27427, prediluted, staining S100A4 in formalin fixed, paraffin embedded melanoma by Immunohistochemistry.



Anti-S100A4 antibody, prediluted (ab27427) at 1/25 dilution + Placenta lysate

Predicted band size: 12 kDa

Observed band size: 12 kDa

Additional bands at: 26 kDa. We are unsure as to the identity of these extra bands.

Western blot - Anti-S100A4 antibody, prediluted (ab27427)

This image refers to an undiluted version of this product.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors