

Product datasheet

Anti-SCGB3A1 antibody ab92350

製品の概要

製品名 Anti-SCGB3A1 antibody  
 製品の詳細 Goat polyclonal to SCGB3A1

**i** This product is a **fast track antibody**. It has been affinity purified and shows high titre values against the immunizing peptide by ELISA.

[Read the terms of use »](#)

種交差性

交差が予測される動物種: Human 

免疫原

Synthetic peptide:  
 NHLIEGSQKCVAE  
 , corresponding to internal sequence amino acids 69 - 81 of Human SCGB3A1 (NP\_443095.2)

 [Run BLAST with](#)

 [Run BLAST with](#)

製品の特性

製品の状態 Liquid  
 保存方法 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.  
 バッファー Preservative: 0.02% Sodium Azide  
 Constituents: 0.5% BSA, Tris buffered saline, pH 7.3  
 精製度 Immunogen affinity purified  
 特記事項 (精製) ab92350 is purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.  
 ポリ/モノ ポリクローナル  
 アイソタイプ IgG

アプリケーション

追加情報 This antibody gave a positive result in ELISA against the immunizing peptide.  
 Not yet tested in other applications.

## ターゲット情報

---

機能	Potential growth inhibitory cytokine.
組織特異性	Highly expressed in breast tissues. Absent in breast cancer cell lines.
配列類似性	Belongs to the secretoglobin family. UGRP subfamily.
細胞内局在	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors