

Product datasheet

Anti-SCGB3A1 antibody ab92350

製品の概要

製品名 Anti-SCGB3A1 antibody

製品の詳細 Goat polyclonal to SCGB3A1

i This product is a **fast track antibody**. It has been affinity purified and shows high titre values against the immunizing peptide by ELISA.

[Read the terms of use »](#)

由来種 Goat

種交差性

交差が予測される動物種: Human

免疫原 Synthetic peptide:
NHLIEGSQKCVAE
, corresponding to internal sequence amino acids 69 - 81 of Human SCGB3A1 (NP_443095.2)

[Run BLAST with](#)

[Run BLAST with](#)

製品の特性

製品の状態 Liquid

保存方法 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

バッファー Preservative: 0.02% Sodium Azide
Constituents: 0.5% BSA, Tris buffered saline, pH 7.3

精製度 Immunogen affinity purified

特記事項 (精製) ab92350 is purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

ポリ/モノ ポリクローナル

アイソタイプ IgG

アプリケーション

Fast track antibodies constitute a diverse group of products that have been released to accelerate your research, but are not yet fully characterized. They have all been affinity purified and show high titre values against the immunizing peptide (by ELISA).

[Fast track terms of use](#)

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB Use at an assay dependent concentration. Predicted molecular weight: 10 kDa.

Preliminary experiments gave an approx 33kDa band in Human Breast cancer and normal Lung lysates after 0.5µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 10.1kDa according to NP_443095.2. The 33kDa band was successfully blocked by incubation with the immunizing peptide.

ターゲット情報

機能	Potential growth inhibitory cytokine.
組織特異性	Highly expressed in breast tissues. Absent in breast cancer cell lines.
配列類似性	Belongs to the secretoglobin family. UGRP subfamily.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors