

Product datasheet

Anti-SLAMF6 antibody [NT-7] ab48444

1 References [画像数 1](#)

製品の概要

| | |
|----------|--|
| 製品名 | Anti-SLAMF6 antibody [NT-7] |
| 製品の詳細 | Mouse monoclonal [NT-7] to SLAMF6 |
| 由来種 | Mouse |
| アプリケーション | 適用あり: Flow Cyt, IP |
| 種交差性 | 交差種: Human |
| 免疫原 | Recombinant fragment from the extracellular domain of Human LY108. |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | Preservative: 0.09% Sodium Azide Constituents: PBS, pH 7.2 |
| 精製度 | Protein G purified |
| ポリ/モノ | モノクローナル |
| クローン名 | NT-7 |
| アイソタイプ | IgG1 |
| 軽鎖の種類 | kappa |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab48444** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| Flow Cyt | | Use 1µg for 10 ⁶ cells. Use 1µg for 10 ⁶ cells in 100 µl volume or 100 µl of whole blood. ab170190 -Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody. |

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| IP | | Use at an assay dependent concentration. |

ターゲット情報

| | |
|-------|--|
| 機能 | Triggers cytolytic activity only in natural killer cells (NK) expressing high surface densities of natural cytotoxicity receptors. |
| 組織特異性 | Expressed by all (resting and activated) natural killer cells (NK), T- and B-lymphocytes. |
| 配列類似性 | Contains 1 Ig-like (immunoglobulin-like) domain. |
| 細胞内局在 | Cell membrane. |

画像



ab48444 staining Human peripheral blood lymphocytes at $1\mu\text{g}$ per 10^6 cells, followed by anti-mouse IgG FITC.

Flow Cytometry - LY108 antibody [NT-7] (ab48444)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors