

## Product datasheet

# Anti-SUC1 antibody [AA10.2] ab24294

### 製品の概要

製品名	Anti-SUC1 antibody [AA10.2]
製品の詳細	Mouse monoclonal [AA10.2] to SUC1
特異性	ab24294 can be used exclusively to capture CDK1 complexes from various cell lines and egg extracts.
アプリケーション	<b>適用あり:</b> WB, IP
種交差性	<b>交差種:</b> Schizosaccharomyces pombe
免疫原	P13 from S.pombe

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	AA10.2
ミエローマ	Sp2
アイソタイプ	IgG1

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab24294** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution. Predicted molecular weight: 13 kDa.
IP		Use at an assay dependent dilution.

## ターゲット情報

---

**関連性** SUC1 is a CDK subunit from *S. pombe* that interacts with CDK1. SUC1 is used exclusively to capture CDK1 complexes from various cell lines and egg extracts. SUC1 belongs to the CKS family. It binds to the catalytic subunit of the cyclin dependent kinase (*cdc2*) and is essential for its biological function. SUC1 forms a homodimer (in presence of zinc).

**細胞内局在** Cytoplasmic and Nuclear

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors