

## Product datasheet

# Anti-Salmonella O3,10 antibody ab80133

### 製品の概要

製品名	Anti-Salmonella O3,10 antibody
製品の詳細	Rabbit polyclonal to Salmonella O3,10
由来種	Rabbit
特異性	Specific to <i>Salmonella</i> O 3,10
アプリケーション	<b>適用あり:</b> Agglutination
種交差性	
免疫原	Tissue/ cell preparation of <i>Salmonella</i> O 3,10 (heat inactivated)
ポジティブ・コントロール	Antigenic suspension
特記事項	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.08% Sodium Azide Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab80133** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Agglutination		Use at an assay dependent concentration. See attached protocol.

## ターゲット情報

---

**関連性** *Salmonella* belongs to the family *Enterobacteriaceae* and is a gram negative non spore forming rod. It causes many types of infections from mild gastroenteritis to life threatening typhoid fever. *Salmonella* can be serotyped according to somatic antigens O (group O: O2-O67 groups), surface Vi, phase I flagellar and phase II flagella and (H types: H-a - H-z81) antigens.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors