

Product datasheet

Anti-Sca1 / Ly6A/E antibody [D7] (PerCP/Cy5.5®) ab93538

画像数 1

製品の概要

製品名	Anti-Sca1 / Ly6A/E antibody [D7] (PerCP/Cy5.5®)
製品の詳細	Rat monoclonal [D7] to Sca1 / Ly6A/E (PerCP/Cy5.5®)
由来種	Rat
標識	PerCP/Cy5.5®. Ex: 482nm, Em: 690nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Mouse
免疫原	IL2 dependent T cell line CTL-L (Mouse).
ポジティブ・コントロール	C57Bl/6 Mouse splenocytes

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 0.1% Gelatin, PBS, 150mM Sodium chloride, pH 7.2
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	D7
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab93538** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		

追加情報

Flow Cyt: Use 0.125 - 0.25 µg for 10<sup>5</sup> - 10<sup>8</sup> cells in a 100 µl volume.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

機能	T-cell activation.
組織特異性	Widely expressed.
配列類似性	Contains 1 UPAR/Ly6 domain.
翻訳後修飾	O-glycosylated. Not N-glycosylated. Not phosphorylated.
細胞内局在	Cell membrane.

## 画像



Flow Cytometry - Anti-Sca1 / Ly6A/E antibody [D7]  
(PerCP/Cy5.5®) (ab93538)

Flow cytometry staining of C57Bl/6 splenocytes with 0.25 µg of PerCP-Cy5.5 Rat IgG2a Isotype Control (open histogram) or 0.125 µg of ab93538 (colored histogram). Cells in the lymphocyte gate were used for analysis.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors