


Product datasheet

Anti-Sin3b antibody ab101841

画像数 2

製品の概要

製品名	Anti-Sin3b antibody
製品の詳細	Rabbit polyclonal to Sin3b
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to a region within amino acids 165-325 of Human Sin3b.
ポジティブ・コントロール	HCT116 whole cell lysate and MCF7 xenograft tissue

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab101841** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 133 kDa.

IHC-P

1/100 - 1/500.

## ターゲット情報


## 関連性

Acts as a transcriptional repressor. Interacts with MYC responsive genes and antagonize MYC oncogenic activities. Interacts with MAD-MAX heterodimers by binding to MAD. The heterodimer then represses transcription by tethering SIN3B to DNA. Also forms a complex with FOXK1 which represses transcription.

## 細胞内局在

Nuclear


## 画像



Western blot - Sin3b antibody (ab101841)

Anti-Sin3b antibody (ab101841) at 1/500 dilution + HCT116 whole cell lysate at 30 µg

**Predicted band size : 133 kDa**



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Sin3b antibody (ab101841)

Immunohistochemical analysis of Sin3b in paraffin embedded MCF7 xenograft, using ab101841 at a 1/100 dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors