

# Anti-Streptococcus Group B II antibody ab78812

## 製品の概要

製品名	Anti-Streptococcus Group B II antibody
製品の詳細	Rabbit polyclonal to Streptococcus Group B II
特異性	ab78812 is specific for Streptococcus Group B II.
アプリケーション	<b>適用あり:</b> Agglutination
種交差性	Reacts with Streptococcus Group B II.
免疫原	Tissue/ cell preparation of Streptococcus Group B II (heat inactivated).
ポジティブ・コントロール	Antigenic suspension.
特記事項	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

## 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.08% Sodium Azide Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab78812** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Agglutination		

## 追加情報

Agg: See attached protocol.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

### 関連性

Hemolytic streptococci are classified into 20 groups, A-V (I and J are omitted), according to serological differences in C-polysaccharide, a specific antigen of the organism. Group B streptococci (GBS), in particular, *Streptococcus agalactiae*, are known to cause serious infectious diseases including neonatal sepsis, pneumonia and meningitis resulting in high fatality rates and sequelae by secondary infection from the mother. GBS serotype is identified by slide agglutination using heat-resistant polysaccharide antigens on the surface of the organism.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors