

Product datasheet

Anti-Streptococcus pneumoniae type 9 antibody ab80486

製品の概要

製品名	Anti-Streptococcus pneumoniae type 9 antibody
製品の詳細	Rabbit polyclonal to Streptococcus pneumoniae type 9
特異性	ab80486 is specific for the capsular antigens of Streptococcus pneumoniae type 9.
アプリケーション	適用あり: Agglutination
種交差性	Reacts with Streptococcus pneumoniae type 9.
免疫原	Tissue/ cell preparation of Streptococcus pneumoniae type 9 (heat inactivated).
ポジティブ・コントロール	Antigenic suspension.
特記事項	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.08% Sodium Azide Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab80486** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Agglutination

追加情報 Agg: See attached protocol.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

関連性

Streptococcus pneumoniae is the most common cause of community-acquired bacterial pneumonia (CAP), and otitis media, and the second most common cause of meningitis in infants and children, which may lead to learning disabilities and hearing impairment. Mortality occurs in ~10-15% of hospitalized adults with invasive disease. While a licensed 23-valent-polysaccharide vaccine is now available, widespread overuse of antibiotics has contributed to the emergence of drug resistant strains. *Streptococcus pneumoniae* is serotyped based on the antigenicity of the capsular polysaccharide.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors