

Product datasheet

Anti-Sua7p antibody ab63909

画像数 1

製品の概要

製品名	Anti-Sua7p antibody
製品の詳細	Rabbit polyclonal to Sua7p
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: <i>Saccharomyces cerevisiae</i>
免疫原	Recombinant full length protein ( <i>Saccharomyces cerevisiae</i> ), purified by chromatography.
ポジティブ・コントロール	<i>Saccharomyces cerevisiae</i> extract.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium Azide Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab63909** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/5000. Detects a band of approximately 41 kDa (predicted molecular weight: 38 kDa).
ELISA		Use at an assay dependent dilution.

## ターゲット情報

**関連性** Sua7p is the TFIIIB of *Saccharomyces cerevisiae*. The fundamental transcription factor TFIIIB has the characteristics of stabilizing the DNA binding of TATA box-binding protein (TBP) and binding directly to DNA by its conformational change. Also its N-terminal region binds to the RNA channel of RNA polymerase undertaking a very important role in the determination of transcription initiation point and promoter clearance.

**細胞内局在** Nuclear

## 画像



Western blot - Anti-Sua7p antibody (ab63909)

**All lanes :** Anti-Sua7p antibody (ab63909) at 1/5000 dilution

**Lane 1 :** Molecular weight markers

**Lane 2 :** *Saccharomyces cerevisiae* extract

**Predicted band size:** 38 kDa

**Observed band size:** 41 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors