

Product datasheet

Anti-Sulphathiazole antibody ab31116

製品の概要

製品名	Anti-Sulphathiazole antibody
製品の詳細	Sheep polyclonal to Sulphathiazole
特異性	% cross reactivity : Sulphathiazole : 100 Sulphapyradine : 0.8 Sulphadiazine : 0.6 Sulphamethizole : <0.2 Sulphachloropyridazine : <0.2 Sulphamerazine : <0.2 Sulphadoxine : <0.2 Sulphaguanidine : <0.2 Sulphisoxazole : <0.2 Sulphanilamide : <0.2 Sulphaquinoxaline : <0.2 Sulphathiazole : <0.2 Sulphamethazine : <0.2 Sulphamethoxypyridazine : <0.2
アプリケーション	適用あり: Competitive ELISA
免疫原	Sulphathiazole conjugated to BSA.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 20mM PBS, 150mM Sodium chloride, pH 7.2
精製度	IgG fraction
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab31116** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Competitive ELISA		1/1000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors