

Product datasheet

Anti-TCHHL1 antibody ab82141

画像数 1

製品の概要

製品名	Anti-TCHHL1 antibody
製品の詳細	Goat polyclonal to TCHHL1
由来種	Goat
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human 交差が予測される動物種: Cow
免疫原	Synthetic peptide: QALEDKQGHPQRE by a Cysteine residue linker, corresponding to internal sequence amino acids 878 - 891 of Human TCHHL1 (NP_001008536.1) Run BLAST with Run BLAST with
ポジティブ・コントロール	Human Skin lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 0.5% BSA, Tris buffered saline, pH 7.3
精製度	Immunogen affinity purified
特記事項 (精製)	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab82141** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション

Abreviews

特記事項

WB

ELISA

追加情報

Peptide ELISA: antibody detection limit dilution 1/4000.


WB: Use at a concentration of 1 - 3 µg/ml. Predicted molecular weight: 99 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

画像



Anti-TCHHL1 antibody (ab82141) at 1 µg/ml +

Human Skin lysate (in RIPA buffer) at 35 µg

Predicted band size: 99 kDa

Observed band size: 99 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors