

## Product datasheet

# Anti-TFIIB antibody [EPR4589] ab109106

リコンビナント RabMAb®

1 References [画像数 1](#)

### 製品の概要

製品名	Anti-TFIIB antibody [EPR4589]
製品の詳細	Rabbit monoclonal [EPR4589] to TFIIB
アプリケーション	<b>適用あり:</b> WB, IP <b>適用なし:</b> Flow Cyt, ICC or IHC-P
種交差性	<b>交差種:</b> Mouse, Rat, Human
免疫原	A synthetic peptide corresponding to the C-terminus of Human TFIIB.
ポジティブ・コントロール	HepG2, HeLa, Jurkat, K562, C6, RAW 264.7, and PC12 cell lysates.
特記事項	This product is a recombinant rabbit monoclonal antibody.  Our RabMAb® technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb® patents</a>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	PBS 49%, Sodium azide 0.01%, Glycerol 50%, BSA 0.05%
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	EPR4589
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab109106** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000. Predicted molecular weight: 35 kDa.
IP		1/10 - 1/100.

**追加情報** Is unsuitable for Flow Cyt, ICC or IHC-P.

### ターゲット情報

<b>機能</b>	General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.
<b>配列類似性</b>	Belongs to the TFIIIB family. Contains 1 TFIIIB-type zinc finger.
<b>細胞内局在</b>	Nucleus.

### 画像



**All lanes** : Anti-TFIIIB antibody [EPR4589]  
(ab109106) at 1/1000 dilution

**Lane 1** : HepG2 cell lysate

**Lane 2** : HeLa cell lysate

**Lane 3** : Jurkat cell lysate

**Lane 4** : K562 cell lysate

**Lane 5** : C6 cell lysate

**Lane 6** : RAW 264.7 cell lysate

**Lane 7** : PC12 cell lysate

Lysates/proteins at 10 µg per lane.

**Predicted band size** : 35 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors