

Product datasheet

Anti-TMEM184B antibody ab106688

画像数 3

製品の概要

製品名	Anti-TMEM184B antibody
製品の詳細	Rabbit polyclonal to TMEM184B
由来種	Rabbit
特異性	ab106688 is predicted to not cross react with other members of the TMEM184 protein family.
アプリケーション	適用あり: WB, IHC-P, ICC/IF
種交差性	交差種: Mouse, Rat, Human
免疫原	A 15 amino acid peptide near the carboxy terminus of Human TMEM184B (NP_036396).
ポジティブ・コントロール	Rat lung tissue lysate; Human lung tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
バッファー	Preservative: 0.02% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項(精製)	Affinity chromatography purified via peptide column.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab106688** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 46 kDa.
IHC-P		Use a concentration of 5 µg/ml.

アプリケーション**Abreviews****特記事項**

ICC/IF

Use a concentration of 20 µg/ml.

ターゲット情報**機能**


May activate the MAP kinase signaling pathway.

配列類似性

Belongs to the TMEM184 family.


細胞内局在

Membrane.

画像

Immunofluorescence of TMEM184B in Human lung tissue with 106688 at 20 µg/mL.

Immunocytochemistry/ Immunofluorescence - Anti-TMEM184B antibody (ab106688)



Lane 1 : Anti-TMEM184B antibody (ab106688) at 1 µg/ml


Lane 2 : Anti-TMEM184B antibody (ab106688) at 2 µg/ml

All lanes : Rat lung tissue lysate

Lysates/proteins at 15 µg per lane.

Western blot - TMEM184B antibody (ab106688)

Predicted band size: 46 kDa



Immunohistochemical analysis of TMEM184B
in Human lung tissue, using ab106688 at 5
 $\mu\text{g/ml}$.

Immunohistochemistry (Formalin/PFA-fixed paraffin-
embedded sections) - TMEM184B antibody
(ab106688)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors