

## Product datasheet

# Anti-TMEM33 antibody ab118435

### 画像数 2

#### 製品の概要

製品名	Anti-TMEM33 antibody
製品の詳細	Rabbit polyclonal to TMEM33
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide conjugated to KLH, from the internal region of Human TMEM33
ポジティブ・コントロール	Human adrenal tissue; Human stomach tissue

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.1% Sodium azide Constituent: 99% PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab118435** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 10 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## ターゲット情報


### 配列類似性

Belongs to the PER33/POM33 family.

### 細胞内局在


Membrane. Melanosome. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

## 画像



ab118435 at 10 µg/ml staining TMEM33 in paraffin embedded formalin fixed Human adrenal tissue by immunohistochemistry (FFPE).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TMEM33 antibody (ab118435)



ab118435 at 10 µg/ml staining TMEM33 in paraffin embedded formalin fixed Human stomach tissue by immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TMEM33 antibody (ab118435)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors