

Product datasheet

Anti-TNFSF18 antibody ab25948

画像数 2

製品の概要

製品名	Anti-TNFSF18 antibody
製品の詳細	Rabbit polyclonal to TNFSF18
由来種	Rabbit
特異性	ab25948 recognises GITRL.
アプリケーション	適用あり: ICC/IF, WB
種交差性	交差種: Mouse, Human
免疫原	Purified recombinant full length protein (Human)
ポジティブ・コントロール	THP-1 cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.02% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25948** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 10 µg/ml.
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 25 kDa (predicted molecular weight: 20 kDa).

ターゲット情報

機能	Cytokine that binds to TNFRSF18/AITR/GITR. Important for interactions between activated T-lymphocytes and endothelial cells and may modulate T-lymphocyte survival in peripheral tissues.
組織特異性	Expressed at high levels in the small intestine, ovary, testis, kidney and endothelial cells.
配列類似性	Belongs to the tumor necrosis factor family.
細胞内局在	Membrane.

画像




All lanes : Anti-TNFSF18 antibody (ab25948)
at 1 μ g/ml

Lane 1 : 5ng of purified recombinant
TNFSF18

Lane 2 : 25ng of purified recombinant
TNFSF18

Lane 3 : 50ng of purified recombinant
TNFSF18

Predicted band size: 20 kDa



ab25948 at 10 μ g/ml staining TNFSF18 in
THP-1 cells by ICC/IF

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors