

Product datasheet

Anti-TRADD antibody [18A11] ab49512

製品の概要

製品名	Anti-TRADD antibody [18A11]
製品の詳細	Mouse monoclonal [18A11] to TRADD
由来種	Mouse
特異性	ab49512 recognises human TNFR1-associated death domain protein (TRADD).
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Prokaryotic recombinant fusion protein corresponding to 85 amino acids of exon 2 of the TRADD.
ポジティブ・コントロール	Immunohistochemistry - EBV-positive Hodgkin's disease.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 15mM Sodium Azide Constituents: Tissue culture supernatant
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	18A11
ミエローマ	P3-NS1/1-Ag4-1
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab49512** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

IHC-P

追加情報	<p>IHC-P: 1/40 - 1/80. High temperature antigen unmasking technique. 60 minutes primary antibody incubation at 25°C.</p> <p>Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.</p>
-------------	--

ターゲット情報	
機能	<p>Adapter molecule for TNFRSF1A/TNFR1 that specifically associates with the cytoplasmic domain of activated TNFRSF1A/TNFR1 mediating its interaction with FADD. Overexpression of TRADD leads to two major TNF-induced responses, apoptosis and activation of NF-kappa-B.</p>
組織特異性	<p>Found in all examined tissues.</p>
配列類似性	<p>Contains 1 death domain.</p>
ドメイン	<p>Requires the intact DEATH domain to associate with TNFRSF1A/TNFR1.</p>
細胞内局在	<p>Cytoplasm > cytoskeleton.</p>

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors