

## Product datasheet

# Anti-TRAM1 antibody ab96106

2 References [画像数 1](#)

### 製品の概要

製品名	Anti-TRAM1 antibody
製品の詳細	Rabbit polyclonal to TRAM1
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide corresponding to a region within amino acids 310-374 (FMMWKFINFQ LRRWREHSAF QAPAVKKKPT VTKGRSSKKG T ENGVNGLT TSNVADSPRN KKEKSS) of Human TRAM1 (NP_055109). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	H1299 whole cell lysate; A431, HeLa and HepG2 lysates

### 法規制情報

**医薬用外毒物**

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab96106** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 43 kDa.

## ターゲット情報

関連性	Function: Stimulatory or required for the translocation of secretory proteins across the ER membrane. Similarity: Belongs to the TRAM family. Contains 1 TLC (TRAM/LAG1/CLN8) domain.
細胞内局在	Endoplasmic reticulum membrane; Multipass membrane protein

## 画像



Anti-TRAM1 antibody (ab96106) at 1/1000 dilution + H1299 cell lysate at 30 µg

**Predicted band size: 43 kDa**

10% SDS-PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors