

Product datasheet

Anti-TRIM58 antibody ab90362

画像数 1

製品の概要

製品名	Anti-TRIM58 antibody
製品の詳細	Rabbit polyclonal to TRIM58
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to a region within internal sequence amino acids 437-486 of Human TRIM58 (NQLFSGLLRP YFFICDATPL ILPPTTIAGS GNWASRDHLD PASDVRDDHL). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	HepG2 cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab90362** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 55 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

配列類似性


Belongs to the TRIM/RBCC family.

Contains 1 B box-type zinc finger.

Contains 1 B30.2/SPRY domain.

Contains 1 RING-type zinc finger.

Anti-TRIM58 antibody 画像



Anti-TRIM58 antibody (ab90362) at 1 μ g/ml +
HepG2 Cell Lysate at 10 μ g

Secondary

HRP conjugated anti-rabbit IgG at 1/50000
dilution

Predicted band size : 55 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors