

Product datasheet

Anti-TRIP antibody ab4533

画像数 2

製品の概要

製品名	Anti-TRIP antibody
製品の詳細	Goat polyclonal to TRIP
由来種	Goat
アプリケーション	適用あり: ICC, WB
種交差性	交差種: Human
免疫原	Synthetic peptide: SLFQAKLDTFLWS , corresponding to C terminal amino acids 457-469 of Human TRIP. Run BLAST with Run BLAST with
特記事項	GenBank Accession Number – NP_005870.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 0.5% BSA, Tris-saline. pH 7.3
精製度	Immunogen affinity purified
特記事項 (精製)	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab4533** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC		Use a concentration of 4 - 6 µg/ml.
WB		Use a concentration of 2 µg/ml. Predicted molecular weight: 53 kDa. additional band of unknown identity was also consistently observed at 38kDa. This band was successfully blocked by incubation with the immunizing peptide.

ターゲット情報

機能	Inhibits activation of NF-kappa-B mediated by TNF.
配列類似性	Contains 1 RING-type zinc finger.
細胞内局在	Cytoplasm. Cytoplasm > perinuclear region.

画像




Anti-TRIP antibody (ab4533) at 2 µg/ml +
Human Testis lysate in RIPA buffer at 35 µg

Developed using the ECL technique.

Predicted band size: 53 kDa

Observed band size: 53 kDa

Additional bands at: 38 kDa. We are unsure as to the identity of these extra bands.



In paraffin embedded Human Spleen shows perinuclear and textured cytoplasmic staining in selective areas.ab4533 (4µg/ml) staining of paraffin embedded Human Spleen. Steamed antigen retrieval with citrate buffer pH 6, HRP-staining.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors