

## Product datasheet

# Anti-Thymosin alpha 1 antibody ab26951

[1 Abreviews](#) [1 References](#)

### 製品の概要

製品名	Anti-Thymosin alpha 1 antibody
製品の詳細	Rabbit polyclonal to Thymosin alpha 1
由来種	Rabbit
特異性	ab26951 recognises Thymosin alpha 1.
アプリケーション	<b>適用あり:</b> ELISA
種交差性	<b>交差種:</b> Rat, Cow, Human
免疫原	Synthetic peptide located within the N terminal half of human Thymosin alpha 1 domain.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab26951** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/100000.

### ターゲット情報

機能	Prothymosin alpha may mediate immune function by conferring resistance to certain opportunistic infections.
配列類似性	Belongs to the pro/parathymosin family.
翻訳後修飾	Covalently linked to a small RNA of about 20 nucleotides.
細胞内局在	Nucleus.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors