

Product datasheet

Anti-Tollip antibody ab62578

★☆☆☆☆ 1 Abreviews 1 References 画像数 2

製品の概要

製品名	Anti-Tollip antibody
製品の詳細	Rabbit polyclonal to Tollip
アプリケーション	適用あり: ELISA, ICC, WB
種交差性	交差種: Mouse, Rat, Human
免疫原	A 16 amino acid peptide from near the center of human Tollip
ポジティブ・コントロール	Rat brain cell lysate and THP-1 cells.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.02% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab62578** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
ICC		Use a concentration of 2 µg/ml.
WB	★☆☆☆☆	Use a concentration of 0.5 - 2 µg/ml. Detects a band of approximately 30 kDa (predicted molecular weight: 30 kDa).

ターゲット情報


機能	Component of the signaling pathway of IL-1 and Toll-like receptors. Inhibits cell activation by microbial products. Recruits IRAK1 to the IL-1 receptor complex. Inhibits IRAK1 phosphorylation and kinase activity.
配列類似性	Belongs to the tollip family. Contains 1 C2 domain. Contains 1 CUE domain.
翻訳後修飾	Phosphorylated by IRAK1 upon stimulation by IL-1 or microbial products.
細胞内局在	Cytoplasm.

画像



Immunocytochemistry - Tollip antibody (ab62578)

ab62578 at 2 µg/ml staining Tollip in THP-1 cells by Immunocytochemistry.



Western blot - Tollip antibody (ab62578)

Lane 1 : Anti-Tollip antibody (ab62578) at 0.5 µg/ml

Lane 2 : Anti-Tollip antibody (ab62578) at 1 µg/ml

Lane 3 : Anti-Tollip antibody (ab62578) at 2 µg/ml

Lane 1 : rat brain cell lysate

Lane 2 : rat brain cell lysate

Lane 3 : rat brain cell lysate

Lysates/proteins at 15 µg per lane.

Predicted band size : 30 kDa

Observed band size : 30 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors